




Saphenous Vein Cannulae

Saphenous Vein Cannulae

▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Ⓢ Saphenous vein cannulae	L201006	25
Outer diameter with the small barb area to the tip: .095 inches ± .006		



Sarns™ INvent Catheters

To drain the left atrium and left ventricle

▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Atrial		
Ⓢ 18 Fr with malleable stylet, 1/4" connector, 18 holes, 14" (35 cm) long	18124	20
18 Fr with soft-flexing stylet, 1/4" connector, 18 holes, 14" (35 cm) long	5859	20
Ventricular		
Ⓢ 18 Fr with malleable stylet, 1/4" connector, 9 holes, 14" (35 cm) long	18080	20
Ⓢ 18 Fr with soft-flexing stylet, 1/4" connector, 9 holes, 14" (35 cm) long	5858	20

✓ Features


Sarns™ INvent Catheters

- Stylet can be shaped
- Choice of malleable or soft-flexing stylet
- Temperature-resistant flexible PVC/polyurethane tubing
- Adjustable stylet length with both extended and recessed positions



Sarns™ Vent Catheters

To vent the left atrium and left ventricle



▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Atrial		
① 18 Fr with 1/4" connector, 18 holes, 14" (35 cm) long	4334	20
Ventricular		
① 18 Fr with 1/4" connector, 9 holes, 14" (35 cm) long	4333	20

✓ Features

Sarns™ Vent Catheters

- Pre-formed removable stylet with retaining latch

Sarns™ Curved Left Heart Vent Catheters

To vent the left atrium and left ventricle


▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
① 20 Fr with 1/4" connector, 18 holes, 13" (33 cm) long	12345	20
① 20 Fr with 1/4" connector, 14 holes, 13" (33 cm) long	14622	20

✓ Features

Sarns™ Curved Left Heart Vent Catheters

- Pre-shaped for easy insertion






1891 Type Left Ventricular Vent Cannulae



1891 Type Left Ventricular Vent Cannulae

▶ Ordering Information

PRESSURE DROP VS. FLOW RATE 1891 TYPE LEFT VENTRICULAR VENT CANNULAE
Terumo in vitro evaluation. Data obtained with water.



DESCRIPTION	CATALOG NO.	UNITS/CASE
22 Fr, 4 side holes, 34 cm (13.4") long ...	L769810

Sarns™ Aortic Needle Vent Catheters

To vent air from the aortic root

▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
17-gauge x 8 mm with 1/4" connector	433620

✓ Features

Sarns™ Aortic Needle Vent Catheters

- Soft, flexible tubing
- Molded tip






1895 Type McKay Left Ventricular Vent Cannulae

1895 Type McKay Left Ventricular Vent Cannulae

▶ Ordering Information

PRESSURE DROP VS. FLOW RATE 1895 TYPE MCKAY LEFT VENTRICULAR VENT CANNULAE


Terumo in vitro evaluation. Data obtained with water.



DESCRIPTION	CATALOG NO.	UNITS/CASE
-------------	-------------	------------

Malleable

18 Fr, 50 cm (20") long	L8581	10
-------------------------	-------	----



Sarns™ Mini Sucker

To clear anastomotic sites of excess blood with gentle suction

▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
-------------	-------------	------------

Mini sucker with 1/4" connector	14991	20
---------------------------------	-------	----

✓ Features

Sarns™ Mini Sucker

- Soft, pliable tip
- Gentle suction