



# Sarns™ Sternal Saw II

A reciprocating saw for rapid and dependable cutting of the sternum

## ▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Sarns™ Sternal Saw II System (220/240 V).....	5590SET	1
includes		
Sternal Saw II (includes saw blade and protector) .....	5590	1
Power unit (motor).....	7084	1
Variable-speed footswitch .....	7085	1
① Flexible drive cable, autoclavable .....	15665	1
<b>Replacement parts for Sarns™ Sternal Saw II System</b>		
① Flexible drive cable .....	15665	1
① Silicone grease for cable .....	131123	1
Cable retainer nut .....	131131	1
<b>Disposables for Sarns™ Sternal Saw II System</b>		
① Sterile blade, non-pyrogenic .....	5589	10
① Non-sterile blade, bulk pack .....	5756	30
Blade protector .....	195135	1
<i>Sarns™ Sternal Saw II Retention System, U.S. patent 5,340,129</i>		
<i>Sarns™ Sternal Saw II Blades, U.S. patent 5,433,457</i>		
<b>Sarns™ Sternal Saw I System</b>		
<b>Replacement parts for Sarns™ Sternal Saw I System</b>		
① Flexible drive cable .....	15665	1
① Silicone grease for cable .....	131123	1
Cable retainer nut .....	131131	1
<b>Disposables for Sarns™ Sternal Saw I System</b>		
① Non-sterile blade with screws .....	6097	6
Blade protector .....	6098	1



Sarns™ Sternal Saw II

## ✓ Features

### Sarns™ Sternal Saw II

- Clear line of sight
- Lightweight design
- Can be ETO sterilized or autoclaved
- Quick-release blades
- Sharp, single-use blades available in sterile or non-sterile bulk pack
- Smooth blade foot piece
- Compatible with original Sarns™ Sternal Saw electro-mechanical drive system